

Adapters for Medical Tap Filter



Type 1 adapter w. stop-valve & 6 L/min. flow restrictor

Item no.: 90-001801
Thread: M24/M22
Material: Brass & WRAS approved plastic



Type 2 adapter w. non-return-valve & 6 L/min. flow restrictor

Item no.: 90-001804
Thread: M24/M22
Material: Brass & WRAS approved plastic



Type 3 adapter w. non-return-valve & stop-valve

Item no.: 90-001805
Thread: M24/M22
Material: Brass & WRAS approved plastic



Type 4 adapter w. non-return-valve, stop-valve & 6 L/min. flow restrictor

Item no.: 90-001806
Thread: M24/M22
Material: Brass & WRAS approved plastic

All adapters are equipped with a M22 adapter ring to ensure installation can be made on either internal or external standard thread.